



TMJ Atlas V.1.0 Installation Instructions for PC and Mac Users

System requirements

512 MB. Ram if Windows XP and 1GB. Ram if Windows Vista

128Mb. graphic card, 3.20 GB. free of hard drive space

1.0 GHz processor or higher

Platform is Windows XP, Vista and Mac OSX

First Steps:

1. Place the Human Occlusion and TMJ DVD into your DVD-Rom

Installation:

1. When the disc is placed into the DVD-Rom reader, one of two things will happen:

PC Users

- a) An installation prompt will appear, and take you through the installation process
- b) If no installation dialogue appears, open "My Computer" then double click on your disc drive and the TMJ Atlas Installer will start

Mac Users

- a) The DVD will mount on your Desktop.
 - b) Double click on the TMJ icon
 - c) Double click on the open box to start TMJ Atlas installation
2. Follow any on screen instructions (you will be asked to enter your system administrator password to continue the installation).
 3. A box will notify you that the installation is complete. Click on "Finish"
 4. To use the TMJ Atlas, click on the Desktop Icon label "Launch TMJ Atlas"